
Evolutis

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**JUST
UNIC®**

Jack-Up **S**houlder implant
for **T**rauma

Evolutis
MOTION INSIDE

Jack-Up Shoulder implant for Trauma

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The JUST UNIC humeral implant is indicated for the repair of 3 and 4 fragments fractures of the proximal humerus independently from the quality of the bone ((1) Report on 26 cases of Bilboquet device at mean FU of 25 months: "the peripheral support on the humeral head optimally counteracts any tendency toward displacement into varus".)

This device allows to realize in a fast and guided way, an anatomical and stable osteosynthesis of the complex fractures of the superior extremity of the humerus ((2) Prospective study of 22 Bilboquet cases at mean FU of 34 months: "Union was obtained in all fractures. There was no secondary tilting of the head, and no migration or pseudarthrosis of the tuberosities").

The impaction of the staple into the cancellous bone of the head provides a firm support which allows for efficient distraction. The stem/sleeve placed in the humeral diaphysis allows for adjustable height distraction of the fracture site to restore an anatomic "Gothic arch" metaphyseal – cephalic position. Subsequent tuberosity re-attachment with sutures is facilitated.

The system was developed by Pr. Levon Doursounian (University Hospital St. Antoine, Paris) in the 1990s allowing for reproducible tuberosity and humeral head consolidation. In cases of symptomatic avascular necrosis of the humeral head, it is possible to protect the osteosynthesis device from the complications related to the necrosis by placing a hollowed head on the Morse taper of the stem.

If the indication justifies the placement of a hollowed head (in particular fracture – dislocations) then the JUST UNIC device allows for the repair in correct position and the healing of the tuberosities.

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(1) A new internal fixation technique for fractures of the proximal humerus - the Bilboquet device: a report on 26 cases.

L.Doursounian et Al. Dept of Orthopaedic Surgery, HEGP, Paris, France. J Shoulder Elbow Surg. July/August 2000: 279-288

(2) Complex proximal humeral fractures: A prospective study of 22 cases treated using the "Bilboquet" device, L. Doursounian, A. Kilinc, B. Cherier, G. Nourissat, Department of Orthopaedic Surgery and Traumatology, Saint-Antoine University Hospital, Paris, France, Orthopaedics & Traumatology: Surgery & Research (2011) 97, 58–66.

- + Strong osteosynthesis
- + Anatomic reduction of the fracture
- + Preservation of the humeral head
- + Healing of the tuberosities
- + Accurate positioning of the height of the humeral head
- + Efficient re-attachment of the tuberosities
- + Auto-graft in contact with the tuberosities (when used with the prosthetic head)

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Humeral head staple

Platform shaped to provide quality support to the femoral head

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Locking screw

Locks the stem to the sleeve at the proper height

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Cementless humeral sleeve

Self-locking in the proximal humeral flare

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Humeral stem

Slides and docks into 1 of 4 positions to the humeral sleeve (3 positions for S6.5S)

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Replacement humeral head

Option when the humeral head has lost its vascularisation. Hollow design for grafting of the humeral head.

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Jack-up forceps: adjust stem level to restore the humeral head height

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